

Peroxidase (POD) Activity Assay Kit (Serum Samples)

Cat No: HR3BC1154

For research use only

Overview

Detection Method	Colorimetric method
Storage	2-8?
Instrument	Spectrophotometer(420 nm)
Assay Time	
Validity	6
Assay Type	Enzyme Activity
Sample Type	Serum
Synonyms	POD
Instrument	Spectrophotometer(420 nm)
Detection Principle	This kit is based on the reaction of hydrogen peroxide catalyzed by peroxidase, it detects the enzymatic activity by measuring the diversification of the absorbency value at 420 nm.
Reagents	
Labware	
Size	50Assays
Sensitivity	0.5?U/mL
Detection Range	0.5-300U/mL
Recovery Rate	