

Glu-Glu Tag (2G1) Mouse mAb

Cat No: HR1AM2508

For research use only

Overview

Product Name	Glu-Glu Tag (2G1) Mouse mAb
Source	Mouse
Applications	WB
Species Reactivity	ALL
Recommended Dilutions	WB 1:2,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1 UCIDU
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	N/A