



ZKSCAN5 Polyclonal Antibody

Cat No: HR1AP6488

For research use only

Overview

Product Name	ZKSCAN5 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	ZKSCAN5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZKSCAN5; KIAA1015; ZFP95; Zinc finger protein with KRAB and SCAN domains 5; Zinc finger protein 95 homolog; Zfp-95
Background	zinc finger with KRAB and SCAN domains 5(ZKSCAN5) Homo sapiens This gene encodes a zinc finger protein of the Kruppel family. The protein contains a SCAN box and a KRAB A domain and may be involved in transcriptional regulation. A similar protein in mouse is differentially expressed in spermatogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015],