

ZN-16 Polyclonal Antibody

Cat No: HR1AP6407

For research use only

Overview

Product Name	ZN-16 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	300kDa
GeneID?Human?	ZNF292
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF292; KIAA0530; Zinc finger protein 292
Background	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 9 C2H2-type zinc fingers.,