

Pan Methylated Lysine Mouse mAb

Cat No: HR1AM2297

For research use only

Overview

Product Name	Pan Methylated Lysine Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	N/A
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	Many proteins are post translationally modified. Modifications such as phosphorylation, glycosylation, ubiquitination and methylation have been shown to play an important role in the development, physiology and disease of animals and plants.