



## TRPS1 Polyclonal Antibody

Cat No: HR1AP8558

For research use only

### Overview

Product Name	TRPS1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	141kDa
GenelD?Human?	TRPS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TRPS1; Zinc finger transcription factor Trps1; Tricho-rhino-phalangeal syndrome type I protein; Zinc finger protein GC79
Background	transcriptional repressor GATA binding 1(TRPS1) Homo sapiens This gene encodes a transcription factor that represses GATA-regulated genes and binds to a dynein light chain protein. Binding of the encoded protein to the dynein light chain protein affects binding to GATA consensus sequences and suppresses its transcriptional activity. Defects in this gene are a cause of tricho-rhino-phalangeal syndrome (TRPS) types I-III. [provided by RefSeq, Jul 2008],