



Treacle Polyclonal Antibody

Cat No: HR1AP8465

For research use only

Overview

Product Name	Treacle Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	152kDa
GenelD?Human?	TCOF1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TCOF1; Treacle protein; Treacher Collins syndrome protein
Background	treacle ribosome biogenesis factor 1(TCOF1) <i>Homo sapiens</i> This gene encodes a nucleolar protein with a LIS1 homology domain. The protein is involved in ribosomal DNA gene transcription through its interaction with upstream binding factor (UBF). Mutations in this gene have been associated with Treacher Collins syndrome, a disorder which includes abnormal craniofacial development. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2008],