



# CMTM1 Polyclonal Antibody

Cat No: HR1AP5752

For research use only

## Overview

Product Name	CMTM1 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
	<a href="http://www.bioelsa.com">www.bioelsa.com</a>
	<a href="mailto:info@bioelsa.com">info@bioelsa.com</a>
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CMTM1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CMTM1; CKLF1; CKLF-like MARVEL transmembrane domain-containing protein 1; Chemokine-like factor superfamily member 1
Background	CKLF like MARVEL transmembrane domain containing 1(CMTM1) Homo sapiens This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and the transmembrane 4 superfamilies of signaling molecules. The protein encoded by this gene may play an important role in testicular development. Alternatively spliced transcript variants encoding different isoforms have been identified. Naturally occurring read-through transcription occurs between this locus and the neighboring locus CKLF (chemokine-like factor).[provided by RefSeq, Feb 2011],