



# IFN- $\alpha$ 1 Polyclonal Antibody

Cat No: HR1AP5011

For research use only

## Overview

[www.bioelsa.com](http://www.bioelsa.com)

[info@bioelsa.com](mailto:info@bioelsa.com)

Product Name	IFN- $\alpha$ 1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	IFNA1/IFNA13
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IFNA1; IFNA13; Interferon alpha-1/13; IFN-alpha-1/13; Interferon alpha-D; LeIF D
Background	interferon alpha 1(IFNA1) Homo sapiens The protein encoded by this gene is produced by macrophages and has antiviral activity. This gene is intronless and the encoded protein is secreted. [provided by RefSeq, Sep 2011],