

# Ephrin-B1/2/3 (phospho Tyr324) Polyclonal Antibody

Cat No: HR1AP6158

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                              |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Ephrin-B1/2/3 (phospho Tyr324) Polyclonal Antibody                                                                                                                                                                                                                                           |
| Source                | Rabbit                                                                                                                                                                                                                                                                                       |
| Applications          | WB,IHC-p,ELISA                                                                                                                                                                                                                                                                               |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                              |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                              |
| Immunogen             |                                                                                                                                                                                                                                                                                              |
| Species               | Rabbit                                                                                                                                                                                                                                                                                       |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                 |
| Isotype               |                                                                                                                                                                                                                                                                                              |
| Clonality             |                                                                                                                                                                                                                                                                                              |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                      |
| Observed band         | 46kDa                                                                                                                                                                                                                                                                                        |
| GeneID?Human?         | EFNB1/EFNB2/EFNB3                                                                                                                                                                                                                                                                            |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                              |
| Cellular localization |                                                                                                                                                                                                                                                                                              |
| Alternative Names     | EFNB1; EFL3; EPLG2; LERK2; Ephrin-B1; EFL-3; ELK ligand; ELK-L; EPH-related receptor tyrosine kinase ligand 2; LERK-2; EFNB2; EPLG5; HTKL; LERK5; Ephrin-B2; EPH-related receptor tyrosine kinase ligand                                                                                     |
| Background            | ephrin B1(EFNB1) Homo sapiens The protein encoded by this gene is a type I membrane protein and a ligand of Eph-related receptor tyrosine kinases. It may play a role in cell adhesion and function in the development or maintenance of the nervous system. [provided by RefSeq, Jul 2008], |