



## TBX22 Polyclonal Antibody

Cat No: HR1AP7601

For research use only

### Overview

Product Name	TBX22 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	58kDa
GenelD?Human?	TBX22
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TBX22; TBOX22; T-box transcription factor TBX22; T-box protein 22
Background	T-box 22(TBX22) <i>Homo sapiens</i> This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental processes. Mutations in this gene have been associated with the inherited X-linked disorder, Cleft palate with ankyloglossia, and it is believed to play a major role in human palatogenesis. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008].