

# SUHW4 Polyclonal Antibody

Cat No: HR1AP7787

For research use only

## Overview

|                       |                                                                                                                         |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name          | SUHW4 Polyclonal Antibody                                                                                               |
| Source                | Rabbit                                                                                                                  |
| Applications          | IHC-p,ELISA                                                                                                             |
| Species Reactivity    | Human                                                                                                                   |
| Recommended Dilutions |                                                                                                                         |
| Immunogen             |                                                                                                                         |
| Species               | Rabbit                                                                                                                  |
| Storage               | -20°C/1 year                                                                                                            |
| Isotype               |                                                                                                                         |
| Clonality             |                                                                                                                         |
| Concentration         | 1 mg/ml                                                                                                                 |
| Observed band         | kDa                                                                                                                     |
| GeneID?Human?         | ZNF280D                                                                                                                 |
| Human Swiss-Prot No.  |                                                                                                                         |
| Cellular localization |                                                                                                                         |
| Alternative Names     | ZNF280D; KIAA1584; SUHW4; ZNF634; Zinc finger protein 280D; Suppressor of hairy wing homolog 4; Zinc finger protein 634 |
| Background            | function:May function as a transcription factor.,similarity:Contains 5 C2H2-type zinc fingers.,                         |