

## TAF I p48 Polyclonal Antibody

Cat No: HR1AP9043

For research use only

### Overview

Product Name	TAF I p48 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53kDa
GeneID?Human?	TAF1A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TAF1A; TATA box-binding protein-associated factor RNA polymerase I subunit A; RNA polymerase I-specific TBP-associated factor 48 kDa; TAFI48; TATA box-binding protein-associated factor 1A; TBP-associa
Background	TATA-box binding protein associated factor, RNA polymerase I subunit A(TAF1A) Homo sapiens This gene encodes a subunit of the RNA polymerase I complex, Selectivity Factor I (SLI). The encoded protein is a TATA box-binding protein-associated factor that plays a role in the assembly of the RNA polymerase I preinitiation complex. Alternate splicing results in multiple transcript variants encoding multiple isoforms.[provided by RefSeq, Jan 2011],