



SMUG1 Polyclonal Antibody

Cat No: HR1AP6475

For research use only

Overview

Product Name	SMUG1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GenID?Human?	SMUG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SMUG1; Single-strand selective monofunctional uracil DNA glycosylase
Background	single-strand-selective monofunctional uracil-DNA glycosylase 1(SMUG1) <i>Homo sapiens</i> This gene encodes a protein that participates in base excision repair by removing uracil from single- and double-stranded DNA. Many alternatively spliced transcript variants exist for this gene; the full-length nature is known for some but not all of the variants. [provided by RefSeq, Aug 2011],