

## Cytochrome C Mouse mAb Cat No: HR1AM2293

For research use only

## Overview

Product Name	Cytochrome C Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-5,000 IHC 1:500-1,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	14kDakDa
GeneID?Human?	54205
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CYC antibody, CYCS antibody, HCS antibody, THC4 antibody, Cytochrome C somatic antibody
Background	The cytochrome complex, or cyt c is a small hemeprotein found loosely associated with the inner membrane of the mitochondrion. It belongs to the cytochrome c family of proteins.