



RNF149 Polyclonal Antibody

Cat No: HR1AP6714

For research use only

Overview

Product Name	RNF149 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	RNF149
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RNF149; DNAPTP2; E3 ubiquitin-protein ligase RNF149; DNA polymerase-transactivated protein 2; RING finger protein 149
Background	domain:The RING-type zinc finger domain mediates binding to an E2 ubiquitin-conjugating enzyme.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 PA (protease associated) domain.,similarity:Contains 1 RING-type zinc finger.,