

PRPF38A Polyclonal Antibody

Cat No: HR1AP8901

For research use only

Overview

Product Name	PRPF38A Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	PRPF38A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PRPF38A; Pre-mRNA-splicing factor 38A
Background	function:May be required for pre-mRNA splicing .,similarity:Belongs to the PRP38 family.,