

PRPF39 Polyclonal Antibody

Cat No: HR1AP7802

For research use only

Overview

Product Name	PRPF39 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	PRPF39
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PRPF39; Pre-mRNA-processing factor 39; PRP39 homolog
Background	function:Involved in pre-mRNA splicing.,similarity:Belongs to the PRP39 family.,similarity:Contains 7 HAT repeats.,