

# POLDIP3 Polyclonal Antibody

Cat No: HR1AP8858

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                  |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | POLDIP3 Polyclonal Antibody                                                                                                                                                                                                                                                                                                                      |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                           |
| Applications          | WB,IHC-p,ELISA                                                                                                                                                                                                                                                                                                                                   |
| Species Reactivity    | Human                                                                                                                                                                                                                                                                                                                                            |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                  |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                  |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                           |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                     |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                  |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                  |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                          |
| Observed band         | 46kDa                                                                                                                                                                                                                                                                                                                                            |
| GeneID?Human?         | POLDIP3                                                                                                                                                                                                                                                                                                                                          |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                  |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                  |
| Alternative Names     | POLDIP3; KIAA1649; PDIP46; Polymerase delta-interacting protein 3; 46 kDa DNA polymerase delta interaction protein; p46; S6K1 Aly/REF-like target; SKAR                                                                                                                                                                                          |
| Background            | DNA polymerase delta interacting protein 3(POLDIP3) Homo sapiens This gene encodes an RRM (RNA recognition motif)-containing protein that participates in the regulation of translation by recruiting ribosomal protein S6 kinase beta-1 to mRNAs. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013], |