

NY-REN-50 Polyclonal Antibody

Cat No: HR1AP7871

For research use only

Overview

Product Name	NY-REN-50 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	ZNF839
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF839; C14orf131; Zinc finger protein 839; Renal carcinoma antigen NY-REN-50
Background	similarity:Contains 1 C2H2-type zinc finger.,