



Ovol1 Polyclonal Antibody

Cat No: HR1AP7562

For research use only

Overview

Product Name	Ovol1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GenelD?Human?	OVOL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	OVOL1; Putative transcription factor Ovo-like 1; hOvo1
Background	ovo like transcriptional repressor 1(OVOL1) Homo sapiens This gene encodes a putative zinc finger containing transcription factor that is highly similar to homologous protein in Drosophila and mouse. Based on known functions in these species, this protein is likely involved in hair formation and spermatogenesis in human as well. [provided by RefSeq, Aug 2011],