



NLE1 Polyclonal Antibody

Cat No: HR1AP7766

For research use only

Overview

Product Name	NLE1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GenID?Human?	NLE1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NLE1; HUSSY-07; Notchless protein homolog 1
Background	similarity:Contains 8 WD repeats.,