




NRIP2 Polyclonal Antibody

Cat No: HR1AP8815

For research use only

Overview

Product Name	NRIP2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	NRIP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NRIP2; Nuclear receptor-interacting protein 2
Background	function:Down-regulates transcriptional activation by nuclear receptors such as NR1F2.,subunit:Interacts with NR1F2, RARA and THRΒ in a ligand-dependent manner.,



bioelsa