

# IP3R-I Polyclonal Antibody

**Cat No: HR1AP9459**

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | IP3R-I Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Applications          | WB,ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Observed band         | 320kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| GeneID?Human?         | ITPR1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Alternative Names     | ITPR1; INSP3R1; Inositol 1,4,5-trisphosphate receptor type 1; IP3 receptor isoform 1; IP3R 1; InsP3R1; Type 1 inositol 1,4,5-trisphosphate receptor; Type 1 InsP3 receptor                                                                                                                                                                                                                                                                                                                                      |
| Background            | inositol 1,4,5-trisphosphate receptor type 1(ITPR1) Homo sapiens This gene encodes an intracellular receptor for inositol 1,4,5-trisphosphate. Upon stimulation by inositol 1,4,5-trisphosphate, this receptor mediates calcium release from the endoplasmic reticulum. Mutations in this gene cause spinocerebellar ataxia type 15, a disease associated with an heterogeneous group of cerebellar disorders. Multiple transcript variants have been identified for this gene. [provided by RefSeq, Nov 2009], |