



Ki-67 Polyclonal Antibody

Cat No: HR1AP3784

For research use only

Overview

Product Name	Ki-67 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	360kDa
GeneID?Human?	MKI67
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MKI67; Antigen KI-67
Background	marker of proliferation Ki-67(MKI67) Homo sapiens This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced transcript variants have been described. A related pseudogene exists on chromosome X. [provided by RefSeq, Mar 2009],