



KKIAMRE Polyclonal Antibody

Cat No: HR1AP9040

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | KKIAMRE Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 56kDa |
| GenID?Human? | CDKL2 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | CDKL2; Cyclin-dependent kinase-like 2; Protein kinase p56 KKIAMRE; Serine/threonine-protein kinase KKIAMRE |
| Background | cyclin dependent kinase like 2(CDKL2) Homo sapiens This gene product is a member of a large family of CDC2-related serine/threonine protein kinases. It accumulates primarily in the cytoplasm, with lower levels in the nucleus. [provided by RefSeq, Jul 2008], |