

KKIAMRE Polyclonal Antibody

Cat No: HR1AP9040

For research use only

Overview

Product Name	KKIAMRE Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	CDKL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDKL2; Cyclin-dependent kinase-like 2; Protein kinase p56 KKIAMRE; Serine/threonine-protein kinase KKIAMRE
Background	cyclin dependent kinase like 2(CDKL2) Homo sapiens This gene product is a member of a large family of CDC2-related serine/threonine protein kinases. It accumulates primarily in the cytoplasm, with lower levels in the nucleus. [provided by RefSeq, Jul 2008],