



HMP19 Polyclonal Antibody

Cat No: HR1AP7660

For research use only

Overview

Product Name	HMP19 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GenelD?Human?	NSG2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NSG2; Neuron-specific protein family member 2; Protein p19; Hmp19
Background	similarity:Belongs to the NSG family.,subcellular location:Mainly Golgi stack, but also found in small vacuolar organelles and multivesicular bodies.,