

## KIAA0100 Polyclonal Antibody

Cat No: HR1AP9199

For research use only

### Overview

Product Name	KIAA0100 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	KIAA0100
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	KIAA0100; BCOX1; UPF0378 protein KIAA0100; Antigen MLLA-22; Breast cancer-overexpressed gene 1 protein
Background	sequence caution:Contaminating sequence.,sequence caution:Sequencing errors.,similarity:Belongs to the UPF0378 family.,tissue specificity:Expressed in pancreas, placenta and up-regulated in breast carcinoma epithelial cells, ductal in situ carcinoma (DCIS), invasive breast carcinoma (IBC) and metastatic breast carcinoma cells (MET).,