

CNTD2 Polyclonal Antibody

Cat No: HR1AP8712

For research use only

Overview

Product Name	CNTD2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CNTD2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CNTD2; Cyclin N-terminal domain-containing protein 2
Background	similarity:Belongs to the cyclin family.,similarity:Contains 1 cyclin N-terminal domain.,