



## Cks2 Polyclonal Antibody

Cat No: HR1AP5771

For research use only

### Overview

Product Name	Cks2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	CKS2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CKS2; Cyclin-dependent kinases regulatory subunit 2; CKS-2
Background	CDC28 protein kinase regulatory subunit 2(CKS2) Homo sapiens CKS2 protein binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function. The CKS2 mRNA is found to be expressed in different patterns through the cell cycle in HeLa cells, which reflects specialized role for the encoded protein. [provided by RefSeq, Jul 2008],