

CDO Polyclonal Antibody

Cat No: HR1AP7600

For research use only

Overview

Product Name	CDO Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CDON
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDON; CDO; Cell adhesion molecule-related/down-regulated by oncogenes
Background	cell adhesion associated, oncogene regulated(CDON) Homo sapiens This gene encodes a cell surface receptor that is a member of the immunoglobulin superfamily. The encoded protein contains three fibronectin type III domains and five immunoglobulin-like C2-type domains. This protein is a member of a cell-surface receptor complex that mediates cell-cell interactions between muscle precursor cells and positively regulates myogenesis. [provided by RefSeq, Aug 2011],