



Cdc5L Polyclonal Antibody

Cat No: HR1AP9229

For research use only

Overview

Product Name	Cdc5L Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	CDC5L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDC5L; KIAA0432; PCDC5RP; Cell division cycle 5-like protein; Cdc5-like protein; Pombe cdc5-related protein
Background	cell division cycle 5 like(CDC5L) Homo sapiens The protein encoded by this gene shares a significant similarity with <i>Schizosaccharomyces pombe</i> cdc5 gene product, which is a cell cycle regulator important for G2/M transition. This protein has been demonstrated to act as a positive regulator of cell cycle G2/M progression. It was also found to be an essential component of a non-snRNA spliceosome, which contains at least five additional protein factors and is required for the second catalytic step of pre-mRNA splicing. [provided by RefSeq, Jul 2008].