




## Cdc50C Polyclonal Antibody

Cat No: HR1AP6740

For research use only

### Overview

Product Name	Cdc50C Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	TMEM30C
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TMEM30C; CDC50C; Cell cycle control protein 50C; Transmembrane protein 30C
Background	TMEM30C, also named as CDC50C, belongs to the CDC50/LEM3 family. TMEM30C is specifically expressed in testis.



**bioelsa**