



## Cdc37L1 Polyclonal Antibody

Cat No: HR1AP7853

For research use only

### Overview

Product Name	Cdc37L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GenID?Human?	CDC37L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDC37L1; CDC37B; HARC; Hsp90 co-chaperone Cdc37-like 1; Hsp90-associating relative of Cdc37
Background	cell division cycle 37 like 1(CDC37L1) Homo sapiens CDC37L1 is a cytoplasmic phosphoprotein that exists in complex with HSP90 (HSPCA; MIM 140571) as well as several other proteins involved in HSP90-mediated protein folding (Scholz et al., 2001 [PubMed 11413142]).[supplied by OMIM, Mar 2008],