



BRWD3 Polyclonal Antibody

Cat No: HR1AP6530

For research use only

Overview

Product Name	BRWD3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	204kDa
GenelD?Human?	BRWD3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BRWD3; Bromodomain and WD repeat-containing protein 3
Background	bromodomain and WD repeat domain containing 3(BRWD3) Homo sapiens The protein encoded by this gene contains a bromodomain and several WD repeats. It is thought to have a chromatin-modifying function, and may thus play a role in transcription. Mutations in this gene cause mental retardation X-linked type 93, which is also referred to as mental retardation X-linked with macrocephaly. This gene is also associated with translocations in patients with B-cell chronic lymphocytic leukemia. [provided by RefSeq, May 2010],