



## MPP9 Polyclonal Antibody

Cat No: HR1AP3927

For research use only

### Overview

Product Name	MPP9 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	78kDa
GenelD?Human?	MPHOSPH9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MPHOSPH9; MPP9; M-phase phosphoprotein 9
Background	PTM:Phosphorylated in M (mitotic) phase.,