



## CdkL4 Polyclonal Antibody

Cat No: HR1AP3056

For research use only

### Overview

Product Name	CdkL4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GenID?Human?	CDKL4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDKL4; Cyclin-dependent kinase-like 4
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,domain:The [NKR]KIAxRE motif seems to be a cyclin-binding region.,similarity:Belongs to the protein kinase superfamily. CMGC Ser/Thr protein kinase family. CDC2/CDKX subfamily.,similarity:Contains 1 protein kinase domain.,