



Cdc50B Polyclonal Antibody

Cat No: HR1AP6077

For research use only

Overview

Product Name	Cdc50B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GenID?Human?	TMEM30B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TMEM30B; CDC50B; Cell cycle control protein 50B; Transmembrane protein 30B
Background	similarity:Belongs to the CDC50/LEM3 family.,

