



BRD3 Polyclonal Antibody

Cat No: HR1AP2901

For research use only

Overview

Product Name	BRD3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	90kDa
GeneID?Human?	BRD3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BRD3; KIAA0043; RING3L; Bromodomain-containing protein 3; RING3-like protein
Background	bromodomain containing 3(BRD3) Homo sapiens This gene was identified based on its homology to the gene encoding the RING3 protein, a serine/threonine kinase. The gene localizes to 9q34, a region which contains several major histocompatibility complex (MHC) genes. The function of the encoded protein is not known. [provided by RefSeq, Jul 2008],