

# Chfr Polyclonal Antibody

Cat No: HR1AP7867

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Chfr Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Applications          | WB,ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Observed band         | 75kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| GeneID?Human?         | CHFR                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Alternative Names     | CHFR; RNF196; E3 ubiquitin-protein ligase CHFR; Checkpoint with forkhead and RING finger domains protein; RING finger protein 196                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Background            | checkpoint with forkhead and ring finger domains(CHFR) Homo sapiens This gene encodes an E3 ubiquitin-protein ligase required for the maintenance of the antephase checkpoint that regulates cell cycle entry into mitosis and, therefore, may play a key role in cell cycle progression and tumorigenesis. The encoded protein has an N-terminal forkhead-associated domain, a central RING-finger domain, and a cysteine-rich C-terminal region. Alternatively spliced transcript variants that encode different protein isoforms have been described. [provided by RefSeq, Mar 2014], |