

NVL Polyclonal Antibody

Cat No: HR1AP4888

For research use only

Overview

Product Name	NVL Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	NVL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NVL; Nuclear valosin-containing protein-like; NVLp; Nuclear VCP-like protein
Background	nuclear VCP-like(NVL) Homo sapiens This gene encodes a member of the AAA (ATPases associated with diverse cellular activities) superfamily. Multiple transcript variants encoding different isoforms have been found for this gene. Two encoded proteins, described as major and minor isoforms, have been localized to distinct regions of the nucleus. The largest encoded protein (major isoform) has been localized to the nucleolus and shown to participate in ribosome biosynthesis (PMID: 15469983, 16782053), while the minor isoform has been localized to the nucleoplasmin. [provided by RefSeq, Aug 2011],