



Lyl-1 Polyclonal Antibody

Cat No: HR1AP7255

For research use only

Overview

Product Name	Lyl-1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenelD?Human?	LYL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LYL1; BHLHA18; Protein lyl-1; Class A basic helix-loop-helix protein 18; bHLHa18; Lymphoblastic leukemia-derived sequence 1
Background	LYL1, basic helix-loop-helix family member(LYL1) Homo sapiens This gene represents a basic helix-loop-helix transcription factor. The encoded protein may play roles in blood vessel maturation and hematopoiesis. A translocation between this locus and the T cell receptor beta locus (GenelD 6957) on chromosome 7 has been associated with acute lymphoblastic leukemia. [provided by RefSeq, Sep 2010],