



## CUL-3 Polyclonal Antibody

Cat No: HR1AP3185

For research use only

### Overview

Product Name	CUL-3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	90kDa
GenID?Human?	CUL3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CUL3; KIAA0617; Cullin-3; CUL-3
Background	cullin 3(CUL3) <i>Homo sapiens</i> This gene encodes a member of the cullin protein family. The encoded protein plays a critical role in the polyubiquitination and subsequent degradation of specific protein substrates as the core component and scaffold protein of an E3 ubiquitin ligase complex. Complexes including the encoded protein may also play a role in late endosome maturation. Mutations in this gene are a cause of type 2E pseudohypoaldosteronism. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Mar 2012].