

Bub3 Polyclonal Antibody

Cat No: HR1AP2911

For research use only

Overview

Product Name	Bub3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	BUB3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BUB3; Mitotic checkpoint protein BUB3
Background	BUB3, mitotic checkpoint protein(BUB3) Homo sapiens This gene encodes a protein involved in spindle checkpoint function. The encoded protein contains four WD repeat domains and has sequence similarity with the yeast BUB3 protein. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008],