



Cyclin C Polyclonal Antibody

Cat No: HR1AP9021

For research use only

Overview

Product Name	Cyclin C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33-37kDa
GenelD?Human?	CCNC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CCNC; Cyclin-C; SRB11 homolog; hSRB11
Background	cyclin C(CCNC) Homo sapiens The protein encoded by this gene is a member of the cyclin family of proteins. The encoded protein interacts with cyclin-dependent kinase 8 and induces the phosphorylation of the carboxy-terminal domain of the large subunit of RNA polymerase II. The level of mRNAs for this gene peaks in the G1 phase of the cell cycle. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008].