

IRF-2 Polyclonal Antibody

Cat No: HR1AP7031

For research use only

Overview

Product Name	IRF-2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	IRF2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IRF2; Interferon regulatory factor 2; IRF-2
Background	interferon regulatory factor 2(IRF2) Homo sapiens IRF2 encodes interferon regulatory factor 2, a member of the interferon regulatory transcription factor (IRF) family. IRF2 competitively inhibits the IRF1-mediated transcriptional activation of interferons alpha and beta, and presumably other genes that employ IRF1 for transcription activation. However, IRF2 also functions as a transcriptional activator of histone H4. [provided by RefSeq, Jul 2008],