



## KKIALRE Polyclonal Antibody

Cat No: HR1AP3788

For research use only

### Overview

Product Name	KKIALRE Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GenelD?Human?	CDKL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDKL1; Cyclin-dependent kinase-like 1; Protein kinase p42 KKIALRE; Serine/threonine-protein kinase KKIALRE
Background	cyclin dependent kinase like 1(CDKL1) Homo sapiens This gene product is a member of a large family of CDC2-related serine/threonine protein kinases. It accumulates primarily in the nucleus. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2013],