



## Cyclin B1 (phospho Ser126) Polyclonal Antibody

Cat No: HR1AP2405

For research use only

### Overview

Product Name	Cyclin B1 (phospho Ser126) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GenelD?Human?	CCNB1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CCNB1; CCNB; G2/mitotic-specific cyclin-B1
Background	cyclin B1(CCNB1) Homo sapiens The protein encoded by this gene is a regulatory protein involved in mitosis. The gene product complexes with p34(cdc2) to form the maturation-promoting factor (MPF). Two alternative transcripts have been found, a constitutively expressed transcript and a cell cycle-regulated transcript, that is expressed predominantly during G2/M phase. The different transcripts result from the use of alternate transcription initiation sites. [provided by RefSeq, Jul 2008],