



## Cdc25C Polyclonal Antibody

Cat No: HR1AP9256

For research use only

### Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Cdc25C Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                            |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                |
| Applications          | IHC-p,IF,ELISA                                                                                                                                                                                                                                                                                                                                                                                                        |
| Species Reactivity    | Human                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                          |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                               |
| Observed band         | kDa                                                                                                                                                                                                                                                                                                                                                                                                                   |
| GeneID?Human?         | CDC25C                                                                                                                                                                                                                                                                                                                                                                                                                |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Alternative Names     | CDC25C; M-phase inducer phosphatase 3; Dual specificity phosphatase Cdc25C                                                                                                                                                                                                                                                                                                                                            |
| Background            | cell division cycle 25C(CDC25C) Homo sapiens This gene encodes a conserved protein that plays a key role in the regulation of cell division. The encoded protein directs dephosphorylation of cyclin B-bound CDC2 and triggers entry into mitosis. It also suppresses p53-induced growth arrest. Multiple alternatively spliced transcript variants of this gene have been described. [provided by RefSeq, Dec 2015], |